



The main goal of Foot Function Lab begins by casting in a non-weight bearing subtalar neutral position. These orthotics are used to balance the musculoskeletal system, control pronation or supination, support the mid metatarsal section, balance the forefoot, and relieve major foot problems such as Heel Spurs, Tendonitis, Fasciitis, etc...

The Custom Easy System is designed to offer a quick diagnosis and present a easy prescription to fit your patient's specific needs and lifestyle. An orthotic may be ordered in any length, form, etc. This allows for ordering the best design and materials to accurately fit and control the patient's foot. System II orthotics are 100% guaranteed to be custom made from each individual's impression, nothing comes off the shelf premade. Orthotics with no heel cups due to precise manufacturing, velcro is included with all models to maintain proper fit without heel cups. These orthotics are manufactured with a metatarsal bar that can be revised to any degree of forefoot control. Any accommodations will be made at no charge for lesions or other conditions that need relief from pain or pressure.

These orthotics are manufactured from a 1-inch and up foam impression box, a plaster cast, or our 3D Light Scanner. Easy System II made with heel cups are made the same as Easy System I.



EZ-1 Standard Spinal Orthotic

Available: Regular or 3/4 length, No Heel Cups, Deluxe Top Cover

Use: Designed for normal daily use, depending on length, in any dress shoe. Used to align the body from heel to toe and to control pronation or supination.

Structure: Manufactured on a flexible thin leather base, cushioned heel pad, with rubber longitudinal arch support, plastazote special beveled metatarsal bar, and standard tan deluxe top cover. And held in place with supplied Velcro® fasteners.

Note: Use regular length for dress shoes, 3/4 length to protect the metatarsal heads when needed, and Full length to protect the toes for sport applications. 3/4 length and Full length extensions have 1/8" Springlite padding under top cover to protect and accommodate metatarsals and toe pressure.



EZ-2 Women's Dress Orthotic SH (Shock Heel)

Available: Regular, 3/4 length, or Full Length, No Heel Cups, Deluxe or 1/16" Pillow Top,

Use: Held with Velcro® fasteners, in all ladies flats, high heels, sling-backs, and sandals. Designed extra narrow to fit the modern contours of ladies dress shoes, is flared out between the 1st and 5th metatarsals for full balance, providing support of the arch, and control pronation or supination of the heel.

Structure: Manufactured on a flexible thin leather base, Deep heel seat of Shock Heel on reinforced base for positive control, with rubber longitudinal arch support, Plastazote special beveled metatarsal bar, and standard deluxe top cover or pillow top. And held in place with supplied Velcro® Fasteners.

Note: A. Regular length is for shoes that fit normally snug leaving no extra room under the metatarsal heads for padding. B. 3/4 Length adds 1/8 of an inch under metatarsal heads requiring shoes to be loose fitting.



EZ-3 Men's Dress Orthotic SH (Shock Heel)

Available: Regular or 3/4 Length, No Heel Cups, Deluxe or 1/16" Pillow Top

Use: Used in all men's dress shoes, loafers, and casual. Designed for the full width contour of men's shoes and is flared out between the 1st and 5th metatarsals for full balance, providing support of the arch, and control pronation or supination of the heel.

Structure: Manufactured on a flexible thin leather base, Deep heel seat of Shock Heel on reinforced base for positive control, with rubber longitudinal arch support, Plastazote special beveled metatarsal bar, and standard deluxe top cover or Pillow Top. And held in place with supplied Velcro® fasteners.

Note: A. Regular length is for shoes that fit normally snug leaving no extra room under the metatarsal heads for padding. B. 3/4 Length adds 1/8 of an inch under metatarsal heads requiring shoes to be loose fitting.