



## SPECIALIZED SPORTS MODELS



### EZ-5A Sports SOCCER

Available: Regular or 3/4 length Deluxe Top Cover or 1/8" Pillow Top

Use: Orthotic has been designed narrow to fit in all soccer shoes, cleats, and running spikes. It has been fitted with a triad metatarsal bar to better performance during activity

Structure: Manufactured on a flexible thin leather base, Deep heel seat of Shock Heel on reinforced base for positive control, with rubber longitudinal arch support, Plastazote special beveled metatarsal bar and standard deluxe top cover or Pillow Top. And held in place with supplied Velcro® fasteners. With No Heel Cups.



### EZ-5B Sports SPIKE

Available: Regular or 3/4 length Deluxe Top Cover or 1/8" Pillow Top

Use: Spike Orthotic has been designed narrow to fit in all soccer shoes, cleats, and running spikes. It has been fitted with a triad metatarsal bar to better performance during activity

Structure: Manufactured on a flexible thin leather base, Deep heel seat of Shock Heel on reinforced base for positive control, with rubber longitudinal arch support, Plastazote special beveled metatarsal bar and standard deluxe top cover or Pillow Top. And held in place with supplied Velcro® fasteners. With No Heel Cups.



### EZ-5C Sports GOLF

Available: 3/4 or Full Length with 1/8" Green Perforated Top

Use: Golf. Extra shock absorption is needed across full bottom of foot surface. Lateral cuboid for better heel control.

Structure: Manufactured on a 1/8" Nickel-Plast base for double shock absorption, SH full heel for extra cushioning, with rubber longitudinal arch support, Plastazote special beveled metatarsal bar, and 1/8" Cushion Green Perforated Cover. Held in place with supplied Velcro® fasteners.

Note: These orthotics must be used in sport shoes with removable inner soles. With No Heel Cups.



### EZ-5D Sports SKATE

Available: Full Length Only

Use: Full-length orthotics designed lightweight. Cushioned orthotic with good flexibility. Built in heel cradle for extra control.

Structure: Manufactured on a 1/8" Nickel-Plast base for double shock absorption, SH full heel for extra cushioning, with rubber longitudinal arch support, Plastazote special beveled metatarsal bar, and 1/8" Cushion Red Perforated Cover. Held in place with supplied Velcro® fasteners.

Note: These orthotics must be used in sport shoes with removable inner soles. With No Heel Cups.