

Custom System Details:

Shell Type

Standard

Low Heel Cup Medial Flange Lateral Flange Medial & Lateral Roberts Plate U.C.B.L. Gait Plate
 Left Right Left Right Left Right Left Right Left Right Left Right Left Right Left Right

Heel Cup Height is approx. 1/4 - 3/8 Heel Cup Height is approx. 1/2 - 5/8 All height and widths are based on each individual foot. Approximations are based on the average foot.

Shell Dimensions

Width:

Narrow - Shell is contained between 1st and 5th met. heads
 Regular - Shell dichotomizes the 1st and 5th met. heads
 Wide - Shell encompasses the 1st and 5th met. heads

Arch Height:

To Cast - Determined during posting to balance the foot
 High - 1/8 of plaster fill
 Medium - 3/8 of plaster fill
 Low - 1/2 of plaster fill

Heel Cup Height:

Low - 3/8
 Medium - 5/8
 High - 3/4

Heel Raise

Left mm Right mm

Pronation Supination

Gait

Shell Padding

Standard, may be added over shell (under vinyl top only)

None
 1/16 Poron®
 1/8 Poron®

Posting Instructions(System I)

Post According To Positive Model(Standard)
 No Posting(neutral shell only)
 Forefoot Only(No measurements needed)
 Post According To Measurements

Left Varus Valgus Right Varus Valgus

RearFoot

ForeFoot

Modifications To Shell

Cut Out 1st Ray In Shell Left Right
 Cut Out 5th Ray In Shell Left Right
 Reinforce Arch Left Right

Modifications To Metatarsal Bar(System II)

Standard: 1/8 Thick 3° Varus Wedge

Changes By Degree: Left Right

Varus

Valgus

Changes By Shape:

Left Right Left Right Left Right

Draw In Any Other Changes:

Left Right

5 4 3 2 1 1 2 3 4 5

Metatarsal Heads Metatarsal Heads

Modifications To Orthotic (System II)

Medial Bars Left Right
 Lateral Bars Left Right

X-13 Heel Lift

Male Female Shoe Size _____

Height: mm

Accommodation Chart:

Mortons Extension

X for R or L Neuroma bar

Mortons Extension

Right Foot **Left Foot**

X for R or L Accommodations

1/8th depth standard or

H Hole

Write in Desired Depth on selected Heads

X for R or L Heel Spur pad

X for Soft U-pad

small medium large
 x for met. pad size and height

small medium large
 x for met. pad size and height

Lateral Spur Pad Medial Spur Pad Center Spur pad Lateral Spur Pad Medial Spur Pad Center Spur pad

Drill Out With Soft Plug Drill Out With Soft Plug