



EZ-4 Western Boot Orthotic SH (Shock Heel)

Available: Regular or 3/4 Length. No Heel Cups, Puff Top Cover

Use: For all western boots. Designed extra narrow between 1st and 5th metatarsal shafts to fit the contour of western boots, it is undercut beneath the arch area, and flared out between 1st and 5th metatarsal heads for full metatarsal control.

Structure: Manufactured on a flexible thin leather base, SH full heel for extra cushioning, with rubber longitudinal arch support, Plastazote special beveled metatarsal bar, and a Microcel Puff top cover for shock absorption and flexibility that won't bottom out. And held in place with supplied Velcro® fasteners.

Note: A. Regular length is for boots that fit normally snug leaving no extra room under the metatarsal heads for padding. B. 3/4 Length adds 1/8 of an inch under metatarsal heads requiring boots to be loose fitting



EZ-5 All Sports Orthotic

Available: 3/4 or Full Length, No Heel Cups, 1/8" Blue Pillow Top

Use: For occupational demanding extra cushioning or normal athletic activities. Golfing, Running, Walking, Soccer, Tennis, Skiing, and Roller Blading.

Structure: Manufactured on a flexible thin leather base, SH full heel for extra cushioning, with rubber longitudinal arch support, Plastazote special beveled metatarsal bar, and Deluxe or Pillow top over 1/16" Springlite heel padding. And held in place with supplied Velcro® fasteners.

Note: A. 3/4 and full length have extra 1/8" springlite padding added to protect metatarsal heads and toes. B. These orthotics must be used in sport shoes with removable inner soles.



EZ-6 High Stress Sports Orthotics

Available: Full Length only, No Heel Cups, with 1/8" Blue Pillow Top

Use: For high impact sports or when extra shock absorption is needed for the full plantar surface of the foot. Such as Aerobics, Basketball, Tennis, Running, and Raquetball.

Structure: Manufactured on a 1/8" Nickel-Plast base for double shock absorption, SH full heel for extra cushioning, with rubber longitudinal arch support, Plastazote special beveled metatarsal bar, and Pillow top cover. And held in place with supplied Velcro® fasteners.

Note: A. These orthotics must be used in sport shoes with removable inner soles.



EZ-7 Cushion Flex Geriatric Orthotic

Available: Full Length Only, No Heel Cups

Use: Full length orthotic designed for the 50+ age group that cannot tolerate a full control sports model. It's a light-weight long lasting cushioned orthotic with good support a flexibility. Good for walking and light sports activities.

Structure: Manufactured on a 1/8" Microcel Puff base for soft support. Beveled Puff heel support to correct pronation or supination, hand-made Puff longitudinal arch support, Plastazote special beveled metatarsal bar, and Puff top cover. And held in place with a supplied Velcro® fasteners.

Note: A. These orthotics must be used in sport shoes with removable inner soles.

X-13 HL Heel Lift



Custom made to any height with Velcro® fasteners.

New Balance Casting Box



Improves foot alignment by up to **80% more!**

NEW 3D SENSE SCANNER!

